








PhaSeal® System Product Overview

The basic PhaSeal product system consists of Protector™, Injector™ and Connector™. A range of products in each category are available to ensure easy and convenient implementation of PhaSeal in the health care setting.



PROTECTOR™

	Denomination / Description	Packaging	
		User box	Transport box
	Protector™ 14 REF P14 Drug vial adapter for closed reconstitution and pressure equalisation. For drug vials with neck Ø13 mm.	50 pcs	200 pcs
	Protector™ 21 REF P21 Drug vial adapter for closed reconstitution and pressure equalisation. For drug vials with neck Ø20 mm.	50 pcs	200 pcs
	Protector™ 28 REF P28 Drug vial adapter for closed reconstitution and pressure equalisation. For drug vials with neck Ø28 mm.	25 pcs	100 pcs
	Protector™ 50 REF P50 Drug vial adapter for closed reconstitution and pressure equalisation. For drug vials with neck Ø20 mm.	25 pcs	100 pcs
	Protector™ 53 REF P53 Drug vial adapter for closed reconstitution and pressure equalisation. For drug vials with neck Ø32 mm.	25 pcs	100 pcs





INJECTOR™

	Denomination / Description	Packaging	
		User box	Transport box
	Injector™ Luer REF N30 Encapsulated cannula that attaches permanently to a disposable syringe with a Luer connection.	50 pcs	200 pcs
	Injector™ Luer Lock REF N35 Encapsulated cannula that attaches to a disposable syringe or IV tubing with a Luer Lock connection.	50 pcs	200 pcs




ADMINISTRATION PRODUCTS – IV LINE ACCESS

	Denomination / Description	Packaging	
		User box	Transport box
	Connector™ Luer Lock REF C40 Attaches to the patient's IV line, providing a closed connection to Injector. Standard Luer Lock connection.	50 pcs	200 pcs
	Y-site Connector™ REF C80 Extension tubing with in-built PhaSeal Connector for closed connection to Injector. No DEHP.	30 pcs	120 pcs

ADMINISTRATION PRODUCTS – IV BAG ACCESS

	Denomination / Description	Packaging	
		User box	Transport box
	Infusion Set REF C50 Standard IV infusion set with in-built PhaSeal Connector for closed addition to IV bag. No DEHP.	–	80 pcs
	Secondary Set REF C61 Secondary set with in-built PhaSeal Connector for closed addition to IV bag. No DEHP.	30 pcs	120 pcs
	Infusion Adapter REF C70 Connection between IV bag and non vented infusion set. In-built PhaSeal Connector for closed addition to IV bag.	–	50 pcs
	Infusion Adapter REF C100 Connection between IV bag and non vented infusion set. In-built PhaSeal Connector for closed addition to IV bag. PVC free.	50 pcs	200 pcs

APPLICATION PRODUCTS

	Denomination / Description	Packaging	
		User box	Transport box
	Assembly Fixture REF M12 Mechanical device that facilitates the attachment of Protector to drug vial. Non sterile.	–	1 pc
	Syringe Tray REF M15 Plastic holder designed to transport filled syringes. Non sterile.	–	200 pcs
	Infusion Clamp REF M25 Clamp for securing Injector-Connector assembly during infusion. Non sterile.	25 pcs	100 pcs